



DDA1 blocking peptide (CDBP5371)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	DDA1 DET1 and DDB1 associated 1 [Homo sapiens (human)]
Official Symbol	DDA1
Synonyms	DDA1; DET1 and DDB1 associated 1; PCIA1; C19orf58; DET1- and DDB1-associated protein 1; PCIA-1; cross-immune reaction antigen PCIA1; placenta cross-immune reaction antigen 1
Entrez Gene ID	79016
mRNA Refseq	NM_024050
Protein Refseq	NP_076955
UniProt ID	Q9BW61