



DDX60 blocking peptide (CDBP5373)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	DDX60 DEAD (Asp-Glu-Ala-Asp) box polypeptide 60 [Homo sapiens (human)]
Official Symbol	DDX60
Synonyms	DDX60; DEAD (Asp-Glu-Ala-Asp) box polypeptide 60; probable ATP-dependent RNA helicase DDX60; DEAD box protein 60
Entrez Gene ID	55601
mRNA Refseq	NM_017631
Protein Refseq	NP_060101
UniProt ID	Q8IY21
Function	ATP binding; double-stranded DNA binding; double-stranded RNA binding; helicase activity;

