



STK17B blocking peptide (CDBP5397)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	STK17B serine/threonine kinase 17b [Homo sapiens (human)]
Official Symbol	STK17B
Synonyms	STK17B; serine/threonine kinase 17b; DRAK2; serine/threonine-protein kinase 17B; death-associated protein kinase-related 2; DAP kinase-related apoptosis-inducing protein kinase 2
Entrez Gene ID	9262
mRNA Refseq	NM_004226
Protein Refseq	NP_004217
UniProt ID	O94768
Function	ATP binding; protein kinase activity; protein serine/threonine kinase activity