



DRGX blocking peptide (CDBP5400)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	DRGX dorsal root ganglia homeobox [Homo sapiens (human)]
Official Symbol	DRGX
Synonyms	DRGX; dorsal root ganglia homeobox; DRG11; PRRXL1; dorsal root ganglia homeobox protein; paired related homeobox-like 1; paired-like homeodomain transcription factor DRG11
Entrez Gene ID	644168
mRNA Refseq	NM_001276451
Protein Refseq	NP_001263380
UniProt ID	A6NNA5
Function	sequence-specific DNA binding; sequence-specific DNA binding transcription factor activity