



CRACR2B blocking peptide (CDBP5404)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	CRACR2B calcium release activated channel regulator 2B [Homo sapiens (human)]
Official Symbol	CRACR2B
Synonyms	CRACR2B; calcium release activated channel regulator 2B; EFCAB4A; EF-hand calcium-binding domain-containing protein 4A; CRAC regulator 2B; CRAC channel regulator 2B; EF-hand calcium binding domain 4A; calcium release-activated calcium channel regulator 2B; Ca2+ release-activated Ca2+ (CRAC) channel regulator 2B
Entrez Gene ID	283229
mRNA Refseq	NM_001286606
Protein Refseq	NP_001273535
UniProt ID	Q8N4Y2

Functioncalcium ion binding
