



CRACR2A blocking peptide (CDBP5405)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	CRACR2A calcium release activated channel regulator 2A [Homo sapiens (human)]
Official Symbol	CRACR2A
Synonyms	CRACR2A; calcium release activated channel regulator 2A; EFCAB4B; EF-hand calcium-binding domain-containing protein 4B; CRAC regulator 2A; CRAC channel regulator 2A; EF-hand calcium binding domain 4B; calcium release-activated calcium channel regulator 2A; Ca ²⁺ release-activated Ca ²⁺ (CRAC) channel regulator 2A
Entrez Gene ID	84766
mRNA Refseq	NM_001144958
Protein Refseq	NP_001138430
UniProt ID	Q9BSW2

Function	calcium ion binding; protein binding
-----------------	--------------------------------------
