



EFHD1 blocking peptide (CDBP5406)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	EFHD1 is an EF-hand domain-containing protein that displays increased expression during neuronal differentiation (Tominaga and Tomooka, 2002 [PubMed 12270117]).[supplied by OMIM, Mar 2008]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	EFHD1 EF-hand domain family, member D1 [Homo sapiens (human)]
Official Symbol	EFHD1
Synonyms	EFHD1; EF-hand domain family, member D1; SWS2; MST133; PP3051; MSTP133; EF-hand domain-containing protein D1; swiprosin-2; EF-hand domain-containing protein 1
Entrez Gene ID	80303
mRNA Refseq	NM_001243252
Protein Refseq	NP_001230181

UniProt ID	Q9BUP0
Function	calcium ion binding