



EFHD2 blocking peptide (CDBP5407)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	EFHD2 EF-hand domain family, member D2 [Homo sapiens (human)]
Official Symbol	EFHD2
Synonyms	EFHD2; EF-hand domain family, member D2; SWS1; EF-hand domain-containing protein D2; swiprosin 1; swiprosin-1; EF hand domain containing 2; EF hand domain family, member D2
Entrez Gene ID	79180
mRNA Refseq	NM_024329
Protein Refseq	NP_077305
UniProt ID	Q96C19
Function	calcium ion binding