



EGI blocking peptide (CDBP5408)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	EGI epilepsy, generalized, idiopathic [Homo sapiens (human)]
Official Symbol	EGI
Synonyms	EGI; epilepsy, generalized, idiopathic; EIG; EIG1
Entrez Gene ID	1957