



KLHL25 blocking peptide (CDBP5414)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	KLHL25 kelch-like family member 25 [Homo sapiens (human)]
Official Symbol	KLHL25
Synonyms	KLHL25; kelch-like family member 25; ENC2; ENC-2; kelch-like protein 25; kelch-like 25; ectodermal-neural cortex 2; BTB/POZ KELCH domain protein; ectoderm-neural cortex protein 2
Entrez Gene ID	64410
mRNA Refseq	NM_022480
Protein Refseq	NP_071925
UniProt ID	Q9H0H3