



ESRP2 blocking peptide (CDBP5422)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ESRP2 is an epithelial cell-type-specific splicing regulator (Warzecha et al., 2009 [PubMed 19285943]).[supplied by OMIM, Aug 2009]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ESRP2 epithelial splicing regulatory protein 2 [Homo sapiens (human)]
Official Symbol	ESRP2
Synonyms	ESRP2; epithelial splicing regulatory protein 2; RBM35B; RNA-binding protein 35B; RNA binding motif protein 35A; RNA binding motif protein 35B; RNA-binding motif protein 35B
Entrez Gene ID	80004
mRNA Refseq	NM_024939
Protein Refseq	NP_079215

UniProt ID	Q9H6T0
Function	mRNA binding; nucleotide binding; poly(A) RNA binding