



FEM1B blocking peptide (CDBP5429)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes an ankyrin repeat protein that belongs to the death receptor-associated family of proteins and plays a role in mediating apoptosis. The encoded protein is also thought to function in the replication stress-induced checkpoint signaling pathway via interaction with checkpoint kinase 1. [provided by RefSeq, Aug 2013]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	FEM1B fem-1 homolog b (C. elegans) [Homo sapiens (human)]
Official Symbol	FEM1B
Synonyms	FEM1B; fem-1 homolog b (C. elegans); F1AA; F1A-ALPHA; FEM1-beta; protein fem-1 homolog B; FEM-1-like death receptor binding protein; fem-1-like in apoptotic pathway protein alpha; fem-1-like death receptor-binding protein alpha
Entrez Gene ID	10116
mRNA Refseq	NM_015322

Protein Refseq	NP_056137
UniProt ID	Q9UK73
Function	death receptor binding; death receptor binding; protein binding; ubiquitin-protein transferase activity