



FAM120A blocking peptide (CDBP5433)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	FAM120A family with sequence similarity 120A [Homo sapiens (human)]
Official Symbol	FAM120A
Synonyms	FAM120A; family with sequence similarity 120A; OSSA; C9orf10; DNATP1; DNATP5; constitutive coactivator of PPAR-gamma-like protein 1; DNA polymerase-transactivated protein 1; DNA polymerase-transactivated protein 5; oxidative stress-associated Src activator; oxidative stress-associated Src activator
Entrez Gene ID	23196
mRNA Refseq	NM_001286722
Protein Refseq	NP_001273651
UniProt ID	Q9NZB2

Functionpoly(A) RNA binding
