



FAM193B blocking peptide (CDBP5434)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	FAM193B family with sequence similarity 193, member B [Homo sapiens (human)]
Official Symbol	FAM193B
Synonyms	FAM193B; family with sequence similarity 193, member B; IRIZIO; protein FAM193B
Entrez Gene ID	54540
mRNA Refseq	NM_001190946
Protein Refseq	NP_001177875
UniProt ID	Q96PV7