



## RMDN2 blocking peptide (CDBP5437)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

### GENE INFORMATION

Gene Name	<a href="#">RMDN2 regulator of microtubule dynamics 2 [ Homo sapiens (human) ]</a>
Official Symbol	RMDN2
Synonyms	RMDN2; regulator of microtubule dynamics 2; RMD2; RMD4; RMD-2; FAM82A; BLOCK18; FAM82A1; PRO34163; PYST9371; regulator of microtubule dynamics protein 2; microtubule-associated protein; family with sequence similarity 82, member A; family with sequence similarity 82, member A1
Entrez Gene ID	<a href="#">151393</a>
mRNA Refseq	<a href="#">NM_001170791</a>
Protein Refseq	<a href="#">NP_001164262</a>
UniProt ID	Q96LZ7