



FCHSD1 blocking peptide (CDBP5441)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	FCHSD1 FCH and double SH3 domains 1 [Homo sapiens (human)]
Official Symbol	FCHSD1
Synonyms	FCHSD1; FCH and double SH3 domains 1; NWK2; FCH and double SH3 domains protein 1; nervous wreck homolog 2
Entrez Gene ID	89848
mRNA Refseq	NM_033449
Protein Refseq	NP_258260
UniProt ID	Q86WN1