



FEZ2 blocking peptide (CDBP5443)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene is an ortholog of the C. elegans unc-76 gene, which is necessary for normal axonal bundling and elongation within axon bundles. Other orthologs include the rat gene that encodes zygin II, which can bind to synaptotagmin. [provided by RefSeq, Jul 2008]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	FEZ2 fasciculation and elongation protein zeta 2 (zygin II) [Homo sapiens (human)]
Official Symbol	FEZ2
Synonyms	FEZ2; fasciculation and elongation protein zeta 2 (zygin II); HUM3CL; fasciculation and elongation protein zeta-2; zygin 2; zygin-2; zygin II; pre-T/NK cell associated protein (3CI)
Entrez Gene ID	9637
mRNA Refseq	NM_001042548
Protein Refseq	NP_001036013

UniProt ID	Q9UHY8
Function	protein binding