



FKBP15 blocking peptide (CDBP5448)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	FKBP15 FK506 binding protein 15, 133kDa [Homo sapiens (human)]
Official Symbol	FKBP15
Synonyms	FKBP15; FK506 binding protein 15, 133kDa; FKBP133; PPP1R76; KIAA0674; FK506-binding protein 15; WAFL; FKBP-15; FKBP-133; 133 kDa FKBP; WASP and FKBP-like protein; FK506-binding protein 133kDa; 133 kDa FK506-binding protein; protein phosphatase 1, regulatory subunit 76
Entrez Gene ID	23307
mRNA Refseq	NM_015258
Protein Refseq	NP_056073
UniProt ID	Q5T1M5

Functionactin binding
