



FNIP2 blocking peptide (CDBP5452)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	FNIP2 folliculin interacting protein 2 [Homo sapiens (human)]
Official Symbol	FNIP2
Synonyms	FNIP2; folliculin interacting protein 2; FNIPL; MAPO1; folliculin-interacting protein 2; FNIP1-like protein; O6-methylguanine-induced apoptosis 1 protein
Entrez Gene ID	57600
mRNA Refseq	NM_020840
Protein Refseq	NP_065891
UniProt ID	Q9P278
Function	protein binding