



FOXRED2 blocking peptide (CDBP5459)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	FOXRED2 FAD-dependent oxidoreductase domain containing 2 [Homo sapiens (human)]
Official Symbol	FOXRED2
Synonyms	FOXRED2; FAD-dependent oxidoreductase domain containing 2; ERFAD; FAD-dependent oxidoreductase domain-containing protein 2; endoplasmic reticulum flavoprotein associated with degradation
Entrez Gene ID	80020
mRNA Refseq	NM_001102371
Protein Refseq	NP_001095841
UniProt ID	Q8IWF2

Function

flavin adenine dinucleotide binding; glycoprotein binding; oxidoreductase activity; protein binding
