



FRMPD1 blocking peptide (CDBP5463)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	FRMPD1 FERM and PDZ domain containing 1 [Homo sapiens (human)]
Official Symbol	FRMPD1
Synonyms	FRMPD1; FERM and PDZ domain containing 1; FRMD2; FERM and PDZ domain-containing protein 1; FERM domain-containing protein 2
Entrez Gene ID	22844
mRNA Refseq	NM_014907
Protein Refseq	NP_055722
UniProt ID	Q5SYB0