



FRMPD3 blocking peptide (CDBP5465)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	FRMPD3 FERM and PDZ domain containing 3 [Homo sapiens (human)]
Official Symbol	FRMPD3
Synonyms	FRMPD3; FERM and PDZ domain containing 3; FERM and PDZ domain-containing protein 3; RP5-1070B1.1
Entrez Gene ID	84443
mRNA Refseq	NM_032428
Protein Refseq	NP_115804
UniProt ID	Q5JV73