



KIAA0195 blocking peptide (CDBP5466)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	KIAA0195 KIAA0195 [Homo sapiens (human)]
Official Symbol	KIAA0195
Synonyms	KIAA0195; TMEM94; uncharacterized protein KIAA0195; transmembrane protein 94
Entrez Gene ID	9772
mRNA Refseq	NM_014738
Protein Refseq	NP_055553
UniProt ID	Q12767