



GMNC blocking peptide (CDBP5477)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	GMNC geminin coiled-coil domain containing [Homo sapiens (human)]
Official Symbol	GMNC
Synonyms	GMNC; geminin coiled-coil domain containing; GEMC1; geminin coiled-coil domain-containing protein 1
Entrez Gene ID	647309
mRNA Refseq	NM_001146686
Protein Refseq	NP_001140158
UniProt ID	A6NCL1
Function	chromatin binding