



Gfra3 blocking peptide (CDBP5480)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	Gfra3 glial cell line derived neurotrophic factor family receptor alpha 3 [Mus musculus (house mouse)]
Official Symbol	Gfra3
Synonyms	Gfra3; glial cell line derived neurotrophic factor family receptor alpha 3; Y15110; GFRalpha3; GDNF family receptor alpha-3; GFR alpha-3; GFR-alpha-3; GDNFR-alpha-3; GDNF receptor alpha-3
Entrez Gene ID	14587
mRNA Refseq	NM_010280
Protein Refseq	NP_034410
UniProt ID	O35118

Function

axon guidance receptor activity; coreceptor activity; protein binding; receptor activity
