



GHDC blocking peptide (CDBP5481)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	GHDC GH3 domain containing [Homo sapiens (human)]
Official Symbol	GHDC
Synonyms	GHDC; GH3 domain containing; LGP1; D11LGP1; GH3 domain-containing protein; homolog of mouse LGP1
Entrez Gene ID	84514
mRNA Refseq	NM_001142622
Protein Refseq	NP_001136094
UniProt ID	Q8N2G8
Function	molecular_function