



GLIPR1L1 blocking peptide (CDBP5486)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	GLIPR1L1 GLI pathogenesis-related 1 like 1 [Homo sapiens (human)]
Official Symbol	GLIPR1L1
Synonyms	GLIPR1L1; GLI pathogenesis-related 1 like 1; PRO7434; ALKN2972; GLIPR1-like protein 1
Entrez Gene ID	256710
mRNA Refseq	NM_152779
Protein Refseq	NP_689992
UniProt ID	Q6UWM5