



GRIK1 blocking peptide (CDBP5501)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Glutamate receptors are the predominant excitatory neurotransmitter receptors in the mammalian brain and are activated in a variety of normal neurophysiologic processes. This gene product belongs to the kainate family of glutamate receptors, which are composed of four subunits and function as ligand-activated ion channels. The subunit encoded by this gene is subject to RNA editing (CAG-> CGG; Q-> R) within the second transmembrane domain, which is thought to alter the properties of ion flow. Alternative splicing, resulting in transcript variants encoding different isoforms, has been noted for this gene. [provided by RefSeq, Jul 2008]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	GRIK1 glutamate receptor, ionotropic, kainate 1 [Homo sapiens (human)]
Official Symbol	GRIK1
Synonyms	GRIK1; glutamate receptor, ionotropic, kainate 1; EAA3; EEA3; GLR5; GLUR5; GluK1; glutamate receptor ionotropic, kainate 1; gluR-5; glutamate receptor 5; excitatory amino acid receptor 3

Entrez Gene ID	2897
mRNA Refseq	NM_000830
Protein Refseq	NP_000821
UniProt ID	P39086
Pathway	Activation of Ca-permeable Kainate Receptor; Activation of Kainate Receptors upon glutamate binding; Activation of Na-permeable Kainate Receptors; Glutamatergic synapse; Ionotropic activity of Kainate Receptors; Neuroactive ligand-receptor interaction; Neuronal System; Neurotransmitter Receptor Binding And Downstream Transmission In The Postsynaptic Cell
Function	extracellular-glutamate-gated ion channel activity; kainate selective glutamate receptor activity