



GRINA blocking peptide (CDBP5508)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	GRINA glutamate receptor, ionotropic, N-methyl D-aspartate-associated protein 1 (glutamate binding) [Homo sapiens (human)]
Official Symbol	GRINA
Synonyms	GRINA; glutamate receptor, ionotropic, N-methyl D-aspartate-associated protein 1 (glutamate binding); LFG1; HNRGW; TMBIM3; NMDARA1; protein lifeguard 1; putative MAPK-activating protein PM02; NMDA receptor glutamate-binding subunit; glutamate [NMDA] receptor-associated protein 1; transmembrane BAX inhibitor motif containing 3; transmembrane BAX inhibitor motif-containing protein 3; glutamate receptor, NMDA subtype, glutamate-binding subunit
Entrez Gene ID	2907
mRNA Refseq	NM_000837

Protein Refseq [NP_000828](#)

UniProt ID Q7Z429
