



HIGD1A blocking peptide (CDBP5531)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	HIGD1A HIG1 hypoxia inducible domain family, member 1A [Homo sapiens (human)]
Official Symbol	HIGD1A
Synonyms	HIGD1A; HIG1 hypoxia inducible domain family, member 1A; HIG1; RCF1a; HIG1 domain family member 1A, mitochondrial; RCF1 homolog A; hypoxia inducible gene 1; HIG1 domain family, member 1A; hypoxia-inducible gene 1 protein
Entrez Gene ID	25994
mRNA Refseq	NM_001099668
Protein Refseq	NP_001093138
UniProt ID	Q9Y241