



HIGD2A blocking peptide (CDBP5532)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	HIGD2A HIG1 hypoxia inducible domain family, member 2A [Homo sapiens (human)]
Official Symbol	HIGD2A
Synonyms	HIGD2A; HIG1 hypoxia inducible domain family, member 2A; RCF1b; HIG1 domain family member 2A, mitochondrial; RCF1 homolog B; HIG1 domain family, member 2A
Entrez Gene ID	192286
mRNA Refseq	NM_138820
Protein Refseq	NP_620175
UniProt ID	Q9BW72