



INSM2 blocking peptide (CDBP5538)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	INSM2 insulinoma-associated 2 [Homo sapiens (human)]
Official Symbol	INSM2
Synonyms	INSM2; insulinoma-associated 2; IA6; IA-6; mlt1; insulinoma-associated protein 2; zinc finger protein IA-6; insulinoma-associated protein IA-6
Entrez Gene ID	84684
mRNA Refseq	NM_032594
Protein Refseq	NP_115983
UniProt ID	Q96T92
Function	DNA binding; metal ion binding