



Dffa blocking peptide (CDBP5540)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	Dffa DNA fragmentation factor, alpha subunit [Mus musculus (house mouse)]
Official Symbol	Dffa
Synonyms	Dffa; DNA fragmentation factor, alpha subunit; ICAD; DFF35; Dff45; ICAD-L; ICAD-S; A330085O09Rik; DNA fragmentation factor subunit alpha; inhibitor of CAD; DNA fragmentation factor 45 kDa subunit
Entrez Gene ID	13347
mRNA Refseq	NM_001025296
Protein Refseq	NP_001020467
UniProt ID	O54786

Pathway	Activation of DNA fragmentation factor; Apoptosis; Apoptosis induced DNA fragmentation; Apoptosis signaling pathway; Apoptotic execution phase; FAS pathway and Stress induction of HSP regulation; PluriNetWork
Function	protein binding