



IFITM1 blocking peptide (CDBP5544)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	IFITM1 interferon induced transmembrane protein 1 [Homo sapiens (human)]
Official Symbol	IFITM1
Synonyms	IFITM1; interferon induced transmembrane protein 1; 9-27; CD225; IFI17; LEU13; DSPA2a; interferon-induced transmembrane protein 1; leu-13 antigen; interferon-induced protein 17; dispanin subfamily A member 2a; interferon-inducible protein 9-27
Entrez Gene ID	8519
mRNA Refseq	NM_003641
Protein Refseq	NP_003632
UniProt ID	P13164

Pathway	Adaptive Immune System; B cell receptor signaling pathway; Cytokine Signaling in Immune system; Immune System; Immunoregulatory interactions between a Lymphoid and a non-Lymphoid cell; Interferon Signaling; Interferon alpha/beta signaling
Function	protein binding; receptor signaling protein activity