



LMNTD1 blocking peptide (CDBP5546)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	LMNTD1 lamin tail domain containing 1 [Homo sapiens (human)]
Official Symbol	LMNTD1
Synonyms	LMNTD1; lamin tail domain containing 1; IFLTD1; PAS1C1; intermediate filament tail domain-containing protein 1; intermediate filament tail domain containing 1
Entrez Gene ID	160492
mRNA Refseq	NM_001145727
Protein Refseq	NP_001139199
UniProt ID	Q8N9Z9
Function	structural molecule activity