



IFNB1 blocking peptide (CDBP5547)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	IFNB1 interferon, beta 1, fibroblast [Homo sapiens (human)]
Official Symbol	IFNB1
Synonyms	IFNB1; interferon, beta 1, fibroblast; IFB; IFF; IFNB; interferon beta; IFN-beta; fibroblast interferon
Entrez Gene ID	3456
mRNA Refseq	NM_002176
Protein Refseq	NP_002167
UniProt ID	P01574
Pathway	Cellular Senescence; Cellular responses to stress; Chagas disease (American

trypanosomiasis); Cytokine Signaling in Immune system; Cytokine-cytokine receptor interaction; Cytokines and Inflammatory Response; Cytosolic DNA-sensing pathway; EBV LMP1 signaling

Function	cytokine activity; cytokine activity; type I interferon receptor binding
-----------------	--
