



IL22 blocking peptide (CDBP5553)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	IL22 interleukin 22 [Homo sapiens (human)]
Official Symbol	IL22
Synonyms	IL22; interleukin 22; TIFa; IL-21; IL-22; ILTIF; IL-TIF; IL-D110; zcyto18; TIFIL-23; interleukin-22; cytokine Zcyto18; IL-10-related T-cell-derived inducible factor; IL-10-related T-cell-derived-inducible factor
Entrez Gene ID	50616
mRNA Refseq	NM_020525
Protein Refseq	NP_065386
UniProt ID	Q9GZX6

Pathway	Cytokine-cytokine receptor interaction; Inflammatory bowel disease (IBD); Jak-STAT signaling pathway
Function	cytokine activity; interleukin-22 receptor binding