



IL31 blocking peptide (CDBP5572)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	IL31 interleukin 31 [Mus musculus (house mouse)]
Official Symbol	IL31
Synonyms	IL31; interleukin 31; 1700013B14Rik; interleukin-31; IL-31
Entrez Gene ID	76399
mRNA Refseq	NM_029594
Protein Refseq	NP_083870
UniProt ID	Q6EAL8
Function	cytokine activity