



BIRC8 blocking peptide (CDBP5575)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	BIRC8 baculoviral IAP repeat containing 8 [Homo sapiens (human)]
Official Symbol	BIRC8
Synonyms	BIRC8; baculoviral IAP repeat containing 8; ILP2; ILP-2; hILP2; baculoviral IAP repeat-containing protein 8; IAP-like protein 2; baculoviral IAP repeat-containing 8; inhibitor of apoptosis-like protein 2; testis-specific inhibitor of apoptosis
Entrez Gene ID	112401
mRNA Refseq	NM_033341
Protein Refseq	NP_203127
UniProt ID	Q96P09

Pathway	Apoptosis; Pathways in cancer; Small cell lung cancer; Toxoplasmosis; Ubiquitin mediated proteolysis
Function	zinc ion binding