



FAM212A blocking peptide (CDBP5582)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	FAM212A family with sequence similarity 212, member A [Homo sapiens (human)]
Official Symbol	FAM212A
Synonyms	FAM212A; family with sequence similarity 212, member A; INKA1; C3orf54; protein FAM212A
Entrez Gene ID	389119
mRNA Refseq	NM_203370
Protein Refseq	NP_976248
UniProt ID	Q96EL1