



ITGA4 blocking peptide (CDBP5584)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The product of this gene belongs to the integrin alpha chain family of proteins. Integrins are heterodimeric integral membrane proteins composed of an alpha chain and a beta chain. This gene encodes an alpha 4 chain. Unlike other integrin alpha chains, alpha 4 neither contains an I-domain, nor undergoes disulfide-linked cleavage. Alpha 4 chain associates with either beta 1 chain or beta 7 chain. [provided by RefSeq, Jul 2008]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ITGA4 integrin, alpha 4 (antigen CD49D, alpha 4 subunit of VLA-4 receptor) [Homo sapiens (human)]
Official Symbol	ITGA4
Synonyms	ITGA4; integrin, alpha 4 (antigen CD49D, alpha 4 subunit of VLA-4 receptor); IA4; CD49D; integrin alpha-4; 269C wild type; integrin alpha 4; integrin alpha-IV; VLA-4 subunit alpha; integrin alpha-4 subunit; CD49 antigen-like family member D; antigen CD49D, alpha-4 subunit of VLA-4 receptor; very late activation protein 4 receptor, alpha 4 subunit

Entrez Gene ID	3676
mRNA Refseq	NM_000885
Protein Refseq	NP_000876
UniProt ID	P13612
Pathway	Adaptive Immune System; Arrhythmogenic right ventricular cardiomyopathy; Arrhythmogenic right ventricular cardiomyopathy (ARVC); Cell adhesion molecules (CAMs); Cell surface interactions at the vascular wall; Dilated cardiomyopathy; ECM-receptor interaction; Extracellular matrix organization
Function	cell adhesion molecule binding; fibronectin binding; metal ion binding; protein binding