



IRF2BP2 blocking peptide (CDBP5600)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes an interferon regulatory factor-2 (IRF2) binding protein that interacts with the C-terminal transcriptional repression domain of IRF2. Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Jul 2008]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	IRF2BP2 interferon regulatory factor 2 binding protein 2 [Homo sapiens (human)]
Official Symbol	IRF2BP2
Synonyms	IRF2BP2; interferon regulatory factor 2 binding protein 2; interferon regulatory factor 2-binding protein 2; IRF-2-binding protein 2
Entrez Gene ID	359948
mRNA Refseq	NM_001077397
Protein Refseq	NP_001070865

UniProt ID	Q7Z5L9
Function	metal ion binding