



IRGC blocking peptide (CDBP5605)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	IRGC immunity-related GTPase family, cinema [Homo sapiens (human)]
Official Symbol	IRGC
Synonyms	IRGC; immunity-related GTPase family, cinema; IFGGE; IRGC1; ligp5; CINEMA; R30953_1; interferon-inducible GTPase 5; interferon inducible GTPase 5; immunity-related GTPase cinema 1; immunity-related GTPase family, cinema 1; interferon-gamma-inducible GTPase IFGGE protein
Entrez Gene ID	56269
mRNA Refseq	NM_019612
Protein Refseq	NP_062558
UniProt ID	Q6NXR0

Function

GTP binding; hydrolase activity, acting on acid anhydrides
