



JAKMIP1 blocking peptide (CDBP5610)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	JAKMIP1 janus kinase and microtubule interacting protein 1 [Homo sapiens (human)]
Official Symbol	JAKMIP1
Synonyms	JAKMIP1; janus kinase and microtubule interacting protein 1; JAMIP1; MARLIN1; Gababrbp; janus kinase and microtubule-interacting protein 1; marlin-1; GABA-B receptor-binding protein; Jak and microtubule interacting protein 1; multiple coiled-coil GABABR1-binding protein; multiple alpha-helices and RNA-linker protein 1
Entrez Gene ID	152789
mRNA Refseq	NM_001099433
Protein Refseq	NP_001092903
UniProt ID	Q96N16

Pathway	Integrated Breast Cancer Pathway
Function	GABA receptor binding; RNA binding; kinase binding; kinesin binding; microtubule binding; protein binding