



JMJD8 blocking peptide (CDBP5622)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	JMJD8 jumonji domain containing 8 [Homo sapiens (human)]
Official Symbol	JMJD8
Synonyms	JMJD8; jumonji domain containing 8; PP14397; C16orf20; jmjC domain-containing protein 8; jumonji domain-containing protein 8
Entrez Gene ID	339123
mRNA Refseq	NM_001005920
Protein Refseq	NP_001005920
UniProt ID	Q96S16