



KANK2 blocking peptide (CDBP5632)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	KANK2 KN motif and ankyrin repeat domains 2 [Homo sapiens (human)]
Official Symbol	KANK2
Synonyms	KANK2; KN motif and ankyrin repeat domains 2; SIP; MXRA3; ANKRD25; KN motif and ankyrin repeat domain-containing protein 2; SRC-interacting protein; SRC1-interacting protein; ankyrin repeat domain 25; SRC-1 interacting protein; SRC-1-interacting protein; matrix-remodelling associated 3; matrix-remodeling-associated protein 3; kidney ankyrin repeat-containing protein 2; ankyrin repeat domain-containing protein 25
Entrez Gene ID	25959
mRNA Refseq	NM_001136191
Protein Refseq	NP_001129663

UniProt ID	Q63ZY3
Function	protein binding