



KANK3 blocking peptide (CDBP5633)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	KANK3 KN motif and ankyrin repeat domains 3 [Homo sapiens (human)]
Official Symbol	KANK3
Synonyms	KANK3; KN motif and ankyrin repeat domains 3; ANKRD47; KN motif and ankyrin repeat domain-containing protein 3; ankyrin repeat domain 47; kidney ankyrin repeat-containing protein 3; ankyrin repeat domain-containing protein 47
Entrez Gene ID	256949
mRNA Refseq	NM_198471
Protein Refseq	NP_940873
UniProt ID	Q6NY19