



NKIRAS2 blocking peptide (CDBP5634)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	NKIRAS2 NFKB inhibitor interacting Ras-like 2 [Homo sapiens (human)]
Official Symbol	NKIRAS2
Synonyms	NKIRAS2; NFKB inhibitor interacting Ras-like 2; KBRAS2; kappaB-Ras2; NF-kappa-B inhibitor-interacting Ras-like protein 2; kappa B-Ras protein 2; I-kappa-B-interacting Ras-like protein 2; NFKB inhibitor interacting Ras-like protein 2
Entrez Gene ID	28511
mRNA Refseq	NM_001001349
Protein Refseq	NP_001001349
UniProt ID	Q9NYR9

Pathway	TNF-alpha/NF-kB Signaling Pathway
Function	GTP binding; GTPase activity