



KCTD15 blocking peptide (CDBP5637)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	KCTD15 potassium channel tetramerization domain containing 15 [Homo sapiens (human)]
Official Symbol	KCTD15
Synonyms	KCTD15; potassium channel tetramerization domain containing 15; BTB/POZ domain-containing protein KCTD15; potassium channel tetramerisation domain containing 15; potassium channel tetramerization domain-containing protein 15
Entrez Gene ID	79047
mRNA Refseq	NM_001129994
Protein Refseq	NP_001123466
UniProt ID	Q96SI1