



KHDC1 blocking peptide (CDBP5639)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	KHDC1 KH homology domain containing 1 [Homo sapiens (human)]
Official Symbol	KHDC1
Synonyms	KHDC1; KH homology domain containing 1; NDG1; C6orf148; bA257K9.4; KH homology domain-containing protein 1
Entrez Gene ID	80759
mRNA Refseq	NM_001251874
Protein Refseq	NP_001238803
UniProt ID	Q4VXA5
Function	RNA binding